



Europa Bioproducts Ltd

15-17 North Street • Wicken
Ely • Cambridgeshire • CB7 5XW
Tel: 0044 (0)1353-721118
Fax: 0044 (0)1353-624589

Bacstain-CTC Rapid Staining Kit

Item#	Unit Size	Highly sensitive fluorescence detection No washing required
BS01-10	100 assays	
BS02-10	100 assays	

Storage Condition: 0-5°C, protect from light

Shipping Condition: ambient temperature

Required Equipment and Materials

10 µl, 1000 µl pipettes, incubator, Microscope (blue excitation filter and red emission filter) or flow cytometer (488 nm blue laser)

Product Description

Colony formation using an agar plate is a very common and reliable method for counting bacterial cells. However, it takes quite a long time to form a colony. Therefore, alternative detection methods have been developed. Bacteria-specific gene amplification methods such as PCR, LAMP, and nucleus staining are quite rapid, but these methods count dead bacteria as well. Therefore, detection of live cell functions is essential to determining the actual number of living bacteria in a sample. Tetrazolium salts can be used to detect respiratory activity of bacterial cells or mitochondria.

CTC is a tetrazolium salt and is reduced by this respiratory activity to form fluorescent CTC formazan on the cell surface. Therefore, CTC is used for specific staining of aerobic live bacteria and can be applied to hard-to-culture bacteria (VNC: viable but non-culturable). CTC forms a fluorescent formazan by an electron transfer system. However, CTC alone is not sensitive enough to stain single cells. Therefore, the CTC-Rapid Staining Kit contains an enhancing reagent that improves the CTC staining efficiency. Compared with staining with CTC only, this staining kit enables rapid and sensitive staining of microorganisms. Maximum wavelengths of the CTC formazan dye are 430 nm or 480 nm for excitation and 630 nm for emission.




Fig. 1 Bacterial cell viability detection mechanism with CTC

Contents of the Kit

CTC Rapid Staining Kit - for Flow cytometry (BS01-10)

CTC...10 mg x 3 vials

Enhancing reagent A ... 100 µl x 1 vial



CERTIFICATE NO. 82422

E-mail: info@europa-bioproducts.com

Web site: <http://www.europa-bioproducts.com>

Registered Office: 15-17 North Street, Wicken, Ely, Cambridgeshire CB7 5XW

Registered Number: 2703381

CTC Rapid Staining Kit- for Microscopy (BS02-10)

CTC...10 mg x 3 vials

Enhancing reagent B ... 500 µl x 1 vial



Europa Bioproducts Ltd

15-17 North Street • Wicken
Ely • Cambridgeshire • CB7 5XW

Tel: 0044 (0)1353-721118
Fax: 0044 (0)1353-624589



CERTIFICATE NO. 32423

E.mail: info@europa-bioproducts.com

Web site: <http://www.europa-bioproducts.com>

Registered Office: 15-17 North Street, Wicken, Ely, Cambridgeshire CB7 5XW

Registered Number: 2703381